


Cardioplegia Cannulae, Retrograde, Manually-Inflating Silicone Balloon, 15 Fr, RCM-14315



The Retrograde Cardioplegia Cannulae are indicated for the infusion of blood or cardioplegia solution into the coronary venous system, intraoperatively. They are available with either a self- or manually-inflating balloon

|                       |           |
|-----------------------|-----------|
| ORDER CODE            | RCM-14315 |
| SIZE                  | 15 Fr     |
| AVAILABLE COATED      | no        |
| EFFECTIVE LENGTH (CM) | 30,5      |
| QTY PER BOX           | 10        |
| LATEX-FREE            | yes       |
| CONNECTOR             | Luer lock |