


Cardioplegia Cannulae, Retrograde, Self-Inflating PVC Balloon, 14 Fr, RCS-13214



The Retrograde Cardioplegia Cannulae are indicated for the infusion of blood or cardioplegia solution into the coronary venous system, intraoperatively. They are available with either a self- or manually-inflating balloon

| | |
|------------------------------|-----------|
| ORDER CODE | RCS-13214 |
| SIZE | 14 Fr |
| AVAILABLE COATED | no |
| EFFECTIVE LENGTH (CM) | 30.5 |
| QTY PER BOX | 5 |
| LATEX-FREE | yes |
| CONNECTOR | Luer lock |