

## MICS Cannulae, Aortic Root Long Needle, 12 Ga (9 Fr), AR-17012



Livanova Aortic Root Long Needle is specifically designed for minimally invasive cardiac surgery in lateral thoracotomy, especially mitral valve procedures. It's longer than a conventional aortic root (35 cm) which allows it to be positioned in lateral thoracotomy.

<b>ORDER CODE</b>	AR-17012
<b>SIZE</b>	12 Ga (9 Fr)
<b>INSERTION NEEDLE</b>	14 Ga - Green
<b>AVAILABLE COATED</b>	no
<b>EFFECTIVE LENGTH (CM)</b>	31,8
<b>QTY PER BOX</b>	10
<b>LATEX-FREE</b>	yes