


# MICS Cannulae, RAP Femoral Venous, 22 Fr Distal/22 Fr Proximal, 200-100



Livanova RAP Femoral Venous Cannulae are designed for improved venous drainage, with a dual stage tip which drains from superior and inferior vena cava. Available in two sizes 22/22 Fr and 23/25 Fr with full-flow drainage capability. A high-quality soft introducer features a smooth transition of the cannula to improve insertion and reduce vascular damage.

Order Code	200-100
Size	22 Fr Distal/22 Fr Proximal
Available Coated	no
Effective Length (cm)	61.5
Qty Per Box	5
Latex-Free	yes