

MICS Cannulae, RAP Femoral Venous, 22 Fr Distal/22 Fr Proximal, 200-100



Livanova RAP Femoral Venous Cannulae are designed for improved venous drainage, with a dual stage tip which drains from superior and inferior vena cava. Available in two sizes 22/22 Fr and 23/25 Fr with full-flow drainage capability. A high-quality soft introducer features a smooth transition of the cannula to improve insertion and reduce vascular damage.

ORDER CODE	200-100
SIZE	22 Fr Distal/22 Fr Proximal
AVAILABLE COATED	no
EFFECTIVE LENGTH (CM)	61,5
QTY PER BOX	5
LATEX-FREE	yes