

# Suction Products, Sumps, Flexible and Malleable Sucker with Fluted Tip, 20 Fr, SU-12502



Suction Wands and Suction Sumps are intended to be used to remove fluids from the surgical field. All of these devices are intended to be used for 6 hours or less. The devices are indicated for use in adult and pediatric patients undergoing surgical procedures requiring cardiopulmonary bypass.

|            |          |
|------------|----------|
| Order Code | SU-12502 |
|------------|----------|

|      |       |
|------|-------|
| Size | 20 Fr |
|------|-------|

|                  |    |
|------------------|----|
| Available Coated | no |
|------------------|----|

|                       |      |
|-----------------------|------|
| Effective Length (cm) | 36.8 |
|-----------------------|------|

|             |    |
|-------------|----|
| Qty Per Box | 10 |
|-------------|----|

|            |     |
|------------|-----|
| Latex-Free | yes |
|------------|-----|