

# Suction Products, Sumps, Flexible and Malleable Sump with Wire-wound Tip, 1/8", SU-20601



Livanova Suckers and Sumps are intended for use to remove excess fluid from the surgical field. They are available with flexible or malleable tubing, with weighted metal spring wound tip that terminates in a rounded metal tip

<b>ORDER CODE</b>	SU-20601
<b>SIZE</b>	1/8 in
<b>AVAILABLE COATED</b>	no
<b>EFFECTIVE LENGTH (CM)</b>	36,8
<b>QTY PER BOX</b>	10
<b>LATEX-FREE</b>	yes