

# Vents, PVC Bullet Tip, without Stylet with Venting Luer Lock, 16 Fr, VT-89616



The Vent Catheters are indicated for use during cardiopulmonary bypass surgery to vent/drain the left ventricle.

<b>ORDER CODE</b>	VT-89616
<b>SIZE</b>	16 Fr
<b>AVAILABLE COATED</b>	no
<b>EFFECTIVE LENGTH (CM)</b>	38
<b>QTY PER BOX</b>	10
<b>LATEX-FREE</b>	yes