

# Vents, PVC Open Blunt Tip, with removeable Stylet, with Two Positions Latching Device, 12 Fr, VT-31212



The Vent Catheters are indicated for use during cardiopulmonary bypass surgery to vent/drain the left ventricle.

<b>ORDER CODE</b>	VT-31212
<b>SIZE</b>	12 Fr
<b>AVAILABLE COATED</b>	no
<b>EFFECTIVE LENGTH (CM)</b>	27,2
<b>QTY PER BOX</b>	10
<b>LATEX-FREE</b>	yes