

Venous Return Cannulae, Single Stage, Right Angle Lighthouse Tip,
Wire-reinforced Tubing, 20 Fr, RV-41020



The Venous Return Cannula is indicated for use in drainage of the superior and inferior vena cava during cardiopulmonary bypass surgery. The Catheter consists of a flexible polyvinylchloride tubing wire reinforced tapered to a flow tip. The Catheters are available with a 90-degree bended distal section with a fluted distal flow tip.

ORDER CODE	RV-41020
SIZE	20 Fr
AVAILABLE COATED	no
QTY PER BOX	10
LATEX-FREE	yes
CONNECTOR ACCEPTANCE	1/4"-3/8