

Vents, PVC Straight, Bullet Tip, with Guidewire Stylet, 13 Fr, VT-84413



The Vent Catheters are indicated for use during cardiopulmonary bypass surgery to vent/drain the left ventricle.

| | |
|------------------------------|----------|
| ORDER CODE | VT-84413 |
| SIZE | 13 Fr |
| AVAILABLE COATED | no |
| EFFECTIVE LENGTH (CM) | 31,8 |
| QTY PER BOX | 10 |
| LATEX-FREE | yes |